



The new Dimension in Video Images


NEW: Sim.CAM 4K Wireless* –

More Details. Better Treatment.

Fascinating images—with clear, sharp contours and strikingly improved color. Details that Full HD was unable to depict are now clearly and distinctly visualized in 4K ULTRA HD. If needed, a cable version of Sim.CAM 4K as a separate arm solution is available.


| Sim.CAM 4K and Sim.CAM 4K Wireless* | |
|-------------------------------------|---|
| Resolution | 8,847,360 pixels |
| 4K Ultra HD | 4096 x 2160 |
| Zoom | 12x optical 12x digital |
| Image positioning | > 360°, mechanically and electronically |
| Video outputs | HDMI 1.4, DVI |

Technical specifications are subject to change
 * Available starting 3rd quarter/2018



Use 4K ULTRA HD in surgical settings for detail-optimized imagery. The high image quality offers you finely structured images with amazing details and thus a brilliant visual representation.



NEW: Sim.CAM HD Wireless – **Full HD 1080** 

More Details. Great Flexibility.

With the new Sim.CAM HD Wireless you obtain a premium HD camera that sends images to your monitor without additional wiring, thanks to the new wireless technology. The integrated HD camera transfers all data sharply and precisely from every selected position. If needed, it is also available as a cable-based solution.

| Sim.CAM HD and Sim.CAM HD Wireless | |
|------------------------------------|---|
| Resolution | 2,000,000 pixels |
| Full HD | 1,920 x 1,080 |
| Aspect ratio | 16:9 |
| Zoom | 10x optical 12x digital |
| Image positioning | > 360°, mechanically and electronically |
| Video outputs | HDMI, DVI |

Technical specifications are subject to change