



Sim.CAM 4K and 4K Wireless\*  
More Details. Better Treatment.



PUSHING TECHNOLOGY TO EXCELLENCE

# The new dimension in video images

Fascinating images – with clear, sharp contours and strikingly improved color. Details that Full HD was unable to depict are now clear and distinct with 4K ULTRA HD.



The Sim.CAM 4K operating panel



External camera

Perfectly integrated, highly flexible, easy-to-clean surfaces




4K technology

Improves the sharpness, color, and contrast of your surgical images

Easy to operate

Camera functions adjustable right on the camera arm



Use 4K in surgical settings for detail-optimized imagery. Its exceptional picture quality provides finely structured images with more accurate details, and thus more brilliant visual representations.

Sim.CAM 4K and 4K Wireless*	
Resolution	approx. 8,510,000 pixels
4K ULTRA HD	3,840 x 2,160
Zoom	20x optical 12x digital
Image reorientation	Limitless, mechanical and electronically
Video outputs	HDMI 1.4 DVI

Tolerance ±10 %; technical specifications are subject to change



PUSHING TECHNOLOGY TO EXCELLENCE