


SIMEON CORONA PACKAGE FOR THE ICU

Easy to install, cost efficient and ready to go – with the SIMEON Corona package for the ICU you can build and expand intensive care capacities to fight against Corona.


ICU Package I – wall version



consisting of:

- Wall mounted ICU pendant
- Sim.LED 250 examination light

ICU Package II – ceiling version



consisting of:

- Ceiling mounted ICU pendant
- Sim.LED 250 examination light

Perform interventions in the ICU to avoid the transport of Corona-patients:

To minimize the infection risk for staff and other patients, Corona patients should be transported as little as possible. To be able to perform small interventions directly in the ICU, the SIMEON Corona package can be upgraded with the Sim.LED 350 examination light.



The basic versions of wall and ceiling installed ICU pendants consist of:

Wet side:

- 12 electric sockets (two circuits)
- 6 earth pin
- 2 RJ45 data sockets
- Arm for infusion rod

Dry side:

- 6 electric sockets (two circuits)
- 2 earth pin
- 2 RJ45 data sockets
- Gas outlets: 2 O₂, 1 AIR, 1 VAC
- 2 shelves
- Monitor arm
- Sim.LED 250 examination light with rail adaptor

As a experienced hospital solution provider, SIMEON is familiar with the requirements of an ICU. The SIMEON Corona package for ICU includes all equipment that is required for the treatment of corona patients. Due to the modular design, tailor-made solutions can be created easily to meet your specific requirements.

Contact us to receive further details about the SIMEON Corona package for the ICU!